

# SmartMicroScope Viewer Online

## User Name and Password Tips

*Web-based software for chromebook, PC and Mac (no install required)*

### User Name:

1. When registering a single user license, the User Name (email) does not have to be a genuine email address.
2. Use something like “SmartMicroScope1@yourschooldistrictdomain.edu” or “Scope1@yourschooldistrictdomain.edu”
3. Put a sticker or masking tape with the user name on each scope. “Scope1” or just “1” to denote which login should be used with that scope.
4. Write the User Name on the inside of the scope box or on the bottom of the gray stand so if it is forgotten or lost it will be easily recovered.

### Password:

1. Everyone has their favorite passwords. **DO NOT USE THAT FOR THIS.**
2. Use a simple password that is easy to remember such as the name you gave during registration. Use “SmartMicroScope1” or “Scope1” as the password.
3. Write the User Name and Password on the inside of the scope box or on the bottom of the gray stand so if it is forgotten or lost it will be easily recovered.

### Single User License

The license is **NOT** linked to a specific scope. You can log in and use any of the SmartMicroScopes, but if you are using it and someone else logs in **USING THE SAME USER NAME AND PASSWORD** it will log you out.

A different User Name and Password **MUST** be used for each simultaneous user.